

Immunofluorescence Staining of Organoids

Materials:

8-Well μ -Slide [80826](#)

Fluoromount-G™ mounting medium [0100-01](#)

Triton™ X-100 [T8787](#)

DAPI [D1306](#)

PBS [10010056](#)

4% PFA [J19943K2](#)

Bovine Serum Albumin (BSA) [9998](#)

Permeabilisation Buffer: 1X PBS / 0.2% Triton™ X-100

Blocking Buffer: 1X PBS / 1% BSA / 0.3% Triton™ X-100

-
1. Culture organoids in 8-well μ -Slides.
 2. Remove culture media and wash organoids with 200 μ L PBS at room temperature for 5 min.
 3. Remove PBS and fix cells with 4 % PFA for 30 min at 4 °C, with (orbital) shaking at around 450 rpm.
 4. Remove PFA (take care not to lose organoids).
 5. Repeat PBS wash three times, 5 min each, with shaking at 45 rpm on rocker.
Note: this step is critical for the removal of Matrigel.
 5. Add 200 μ L Permeabilisation Buffer and incubate for 30 min at room temperature.
 6. Add 200 μ L Blocking Buffer and incubate for 30 min at room temperature. Prepare primary antibodies during the incubation.
 7. Add 200 μ L primary antibody diluted in Blocking Buffer, and incubate for 2 hours at room temperature or overnight at 4 °C.
 8. Repeat PBS wash three times. Prepare secondary antibodies during the washes.
 9. Add 200 μ L secondary antibodies and DAPI diluted in Blocking Buffer, leave for 1 hr at room temperature.
Note: shield cells from light using foil from now on.
 10. Repeat PBS wash three times.
 11. Wash the cells in ddH₂O.
 12. Mount the cells in Fluoromount-G™ mounting medium, leave at room temperature until dry, and store slide at 4 °C prior to imaging.